

Boston Scoliosis Brace Customised

The Original Boston Brace system for the non-operative treatment of Idiopathic Scoliosis, utilises a prefabricated, standardised symmetrical module chosen on the basis of the patients physical dimensions. The Boston Brace system is still the preferred choice in the non-operative treatment of scoliosis. The brace is manufactured from Orthopaedic grade copolymer with a Hypo-allergenic Polyethylene foam liner.



Recommended Range Of Application

Idiopathic Scoliosis

Measurements

Can be made to measurements or cast. Please use our universal measurement form for customised brace ordering which can be found in our download centre. When ordering Milwaukee Superstructure please contact us for guidance on measurements.

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.

Item No.	Model	Lordosis degree
BB02	Made to Measure Module	Customised
BB03	Made to Cast Module	Customised
BB05	Finishing to First Fitting	Customised
BB06	Milwaukee Superstructure	Customised
BB09	All Over Transfer Pattern	Customised