

# Boston Scoliosis Brace Stock Size 15°

The Original Boston Brace system for the non-operative treatment of Idiopathic Scoliosis, utilises a prefabricated, standardised symmetrical module chosen on the basis of the patients physical dimensions. The Boston Brace system is still the preferred choice in the non-operative treatment of scoliosis. The brace is manufactured from Orthopaedic grade copolymer with a Hypo-allergenic Polyethylene foam liner.

## **Recommended Range Of Application**

Idiopathic Scoliosis

## **Measurements**

Available in 30 stock sizes with either a 0°, 15° or 30° Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree
15J1	Stock Size Module	J1	54	40	15°
15J2	Stock Size Module	J2	58	42	15°
15J3	Stock Size Module	J3	62	45	15°
15J4	Stock Size Module	J4	64	47	15°
15Y1A	Stock Size Module	Y1A	66	51	15°
15Y1	Stock Size Module	Y1	67	46	15°
15Y2	Stock Size Module	Y2	71	47	15°
15Y2A	Stock Size Module	Y2A	72	55	15°
15Y3	Stock Size Module	Y3	75	48	15°
15Y3A	Stock Size Module	Y3A	78	57	15°
15Y4	Stock Size Module	Y4	79	51	15°
15Y4F	Stock Size Module	Y4F	80	57	15°
15Y5	Stock Size Module	Y5	82	52	15°
15Y5F	Stock Size Module	Y5F	82	58	15°
15Y6	Stock Size Module	Y5	84	56	15°
15Y6A	Stock Size Module	Y6A	84	63	15°
15Y7	Stock Size Module	Y7	86	54	15°
15Y7F	Stock Size Module	Y7F	86	60	15°
15Y7A	Stock Size Module	Y7A	86	69	15°
15Y8	Stock Size Module	Y8	88	58	15°
15Y8F	Stock Size Module	Y8F	88	64	15°
15Y8A	Stock Size Module	Y8A	90	72	15°
15Y9	Stock Size Module	Y9	92	60	15°
15Y9A	Stock Size Module	Y9A	92	63	15°
15Y9F	Stock Size Module	Y9F	92	68	15°
15A1	Stock Size Module	A1	94	62	15°
15A2	Stock Size Module	A2	94	70	15°
15A2A	Stock Size Module	A2A	96	66	15°
15A3	Stock Size Module	A3	97	68	15°
15A4	Stock Size Module	A4	97	76	15°