

Boston Soft AFO

The Boston Soft AFO is a resting splint which is an inexpensive alternative to a custom AFO. It is lightweight and manufactured from 3 mm soft-feel polyethylene foam with a 3 mm polypropylene heel insert. This AFO if fully adjustable and can be easily trimmed with scissors.

Recommended Range Of Application

Myelomeningocele, Muscular Dystrophy, Cerebral Palsy, Trauma and Heel Ulceration.

Contraindications

None known.

Measurements

Can be used bilaterally.

Measurement 1: UK Shoe Size.

Measurement 2: European Shoe Size.



Item No.	Model	Size	Measurement 1	Measurement 2
AFO-A	Adult without stabilizer	A	1-2	32-34
AFO-AS	Adult with stabilizer	AS	1-2	32-34
AFO-B	Adult without stabilizer	B	3-4	35-36
AFO-BS	Adult with stabilizer	BS	3-4	35-36
AFO-C	Adult without stabilizer	C	5-6	37-39
AFO-CS	Adult with stabilizer	CS	5-6	37-39
AFO-D	Adult without stabilizer	D	7-8	40-41
AFO-DS	Adult with stabilizer	DS	7-8	40-41
AFO-E	Adult without stabilizer	E	9-10	42-44
AFO-ES	Adult with stabilizer	ES	9-10	42-44
AFO-F	Adult without stabilizer	F	11-12	45-46
AFO-FS	Adult with stabilizer	FS	11-12	45-46
AFO-P1	Paediatric without stabilizer	P1	2-3	19-20
AFO-P1S	Paediatric with stabilizer	P1S	2-3	19-20
AFO-P2	Paediatric without stabilizer	P2	4-5	21-22
AFO-P2S	Paediatric with stabilizer	P2S	4-5	21-22
AFO-P3	Paediatric without stabilizer	P3	6-7	23-24
AFO-P3S	Paediatric with stabilizer	P3S	6-7	23-24
AFO-P4	Paediatric without stabilizer	P4	8-9	25-26
AFO-P4S	Paediatric with stabilizer	P4S	8-9	25-26
AFO-P5	Paediatric without stabilizer	P5	10-11	27-29
AFO-P5S	Paediatric with stabilizer	P5S	10-11	27-29
AFO-P6	Paediatric without stabilizer	P6	12-13	30-31
AFO-P6S	Paediatric with stabilizer	P6S	12-13	30-31
AFO-SK	Stabilizer Kit	-	-	-