

# Boston Soft AFO

The Boston Soft AFO is a resting splint which is an inexpensive alternative to a custom AFO. It is lightweight and manufactured from 3 mm soft-feel polyethylene foam with a 3 mm polypropylene heel insert. This AFO is fully adjustable and can be easily trimmed with scissors.

## Recommended Range Of Application

Myelomeningocele, Muscular Dystrophy, Cerebral Palsy, Trauma and Heel Ulceration.

## Contraindications

None known.

## Measurements

Can be used bilaterally.

Measurement 1: UK Shoe Size.

Measurement 2: European Shoe Size.



Item No.	Model	Size	Measurement 1	Measurement 2
AFO-A	Adult without stabilizer	A	1-2	32-34
AFO-AS	Adult with stabilizer	AS	1-2	32-34
AFO-B	Adult without stabilizer	B	3-4	35-36
AFO-BS	Adult with stabilizer	BS	3-4	35-36
AFO-C	Adult without stabilizer	C	5-6	37-39
AFO-CS	Adult with stabilizer	CS	5-6	37-39
AFO-D	Adult without stabilizer	D	7-8	40-41
AFO-DS	Adult with stabilizer	DS	7-8	40-41
AFO-E	Adult without stabilizer	E	9-10	42-44
AFO-ES	Adult with stabilizer	ES	9-10	42-44
AFO-F	Adult without stabilizer	F	11-12	45-46
AFO-FS	Adult with stabilizer	FS	11-12	45-46
AFO-P1	Paediatric without stabilizer	P1	2-3	19-20
AFO-P1S	Paediatric with stabilizer	P1S	2-3	19-20
AFO-P2	Paediatric without stabilizer	P2	4-5	21-22
AFO-P2S	Paediatric with stabilizer	P2S	4-5	21-22
AFO-P3	Paediatric without stabilizer	P3	6-7	23-24
AFO-P3S	Paediatric with stabilizer	P3S	6-7	23-24
AFO-P4	Paediatric without stabilizer	P4	8-9	25-26
AFO-P4S	Paediatric with stabilizer	P4S	8-9	25-26
AFO-P5	Paediatric without stabilizer	P5	10-11	27-29
AFO-P5S	Paediatric with stabilizer	P5S	10-11	27-29
AFO-P6	Paediatric without stabilizer	P6	12-13	30-31
AFO-P6S	Paediatric with stabilizer	P6S	12-13	30-31
AFO-SK	Stabilizer Kit	-	-	-