

Boston Hip Abduction Brace

The Boston Hip Abduction Brace is a lightweight paediatric orthosis with a frame that is manufactured from sturdy polypropylene and a closed cell, hypoallergenic foam liner which does not retain moisture. This orthosis maintains concentric reduction of hips in patients with dysplasia. The product can replace spica casts and harnesses for certain developmental dysplasia patients.



Recommended Range Of Application

Congenital Hip Dislocation, Hip Dysplasia and Cerebral Palsy.

Contraindications

None known.

Item No.	Model	Size	Pelvic Circumference
HPXXS	Hip Abduction Brace	XX-Small	26-30cm
HPXS	Hip Abduction Brace	X-Small	30-38cm
HPS	Hip Abduction Brace	Small	38-46cm
HPM	Hip Abduction Brace	Medium	46-51cm
HPL	Hip Abduction Brace	Large	51-56cm