

Boston Overlap Brace (BOB) Stock Size 0°

The Original Boston Overlap Brace (BOB) is an anterior opening LSO thermoplastic brace specifically designed for patients suffering from chronic or temporary low back pain.

Recommended Range Of Application

Mechanical back pain, Discongenic low back pain, Compression fractures, Sports injuries, Spondylolysis, Spondylolisthesis, Osteoporosis, Spinal Arthrosis, Spinal Stenosis, Pre and Post-operative care.

Measurements

Available in 20 stock sizes with either 0°, 15° or 30° Lordosis. Measurements in centimetres: Measurement 1: Hip, Measurement 2: Waist, Measurement 3: Chest. This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree
001	Stock Size Module, Unlined	1	71	51	0°
001D	Stock Size Module, Unlined	1D	71	45	0°
002	Stock Size Module, Unlined	2	76	61	0°
002D	Stock Size Module, Unlined	2D	76	51	0°
003	Stock Size Module, Unlined	3	81	71	0°
003D	Stock Size Module, Unlined	3D	81	66	0°
004	Stock Size Module, Unlined	4	86.5	76	0°
004D	Stock Size Module, Unlined	4D	86.5	58.5	0°
005	Stock Size Module, Unlined	5	91.5	81	0°
005D	Stock Size Module, Unlined	5D	91.5	66	0°
006	Stock Size Module, Unlined	6	96.5	86.5	0°
006D	Stock Size Module, Unlined	6D	96.5	74	0°
007	Stock Size Module, Unlined	7	101.5	91.5	0°
007D	Stock Size Module, Unlined	7D	101.5	81	0°
008	Stock Size Module, Unlined	8	109	96.5	0°
008D	Stock Size Module, Unlined	8D	109	84	0°
009	Stock Size Module, Unlined	9	114	100	0°
009D	Stock Size Module, Unlined	9D	114	87	0°
0010	Stock Size Module, Unlined	10	117	104	0°
0010D	Stock Size Module, Unlined	10D	117	96	0°
001L	Stock Size Module, Lined	1	71	51	0°
001DL	Stock Size Module, Lined	1D	71	45	0°
002L	Stock Size Module, Lined	2	76	45	0°
002DL	Stock Size Module, Lined	2D	76	51	0°
003L	Stock Size Module, Lined	3	81	71	0°
003DL	Stock Size Module, Lined	3D	81	66	0°
004L	Stock Size Module, Lined	4	86.5	76	0°
004DL	Stock Size Module, Lined	4D	86.5	58.5	0°
005L	Stock Size Module, Lined	5	91.5	81	0°
005DL	Stock Size Module, Lined	5D	91.5	66	0°
006L	Stock Size Module, Lined	6	96.5	86.5	0°
006DL	Stock Size Module, Lined	6D	96.5	74	0°
007L	Stock Size Module, Lined	7	101.5	91.5	0°
007DL	Stock Size Module, Lined	7D	101.5	81	0°
008L	Stock Size Module, Lined	8	109	96.5	0°
008DL	Stock Size Module, Lined	8D	109	84	0°
009L	Stock Size Module, Lined	9	114	100	0°
009DL	Stock Size Module, Lined	9D	114	87	0°
0010L	Stock Size Module, Lined	10	117	104	0°
0010DL	Stock Size Module, Lined	10D	117	96	0°