

Boston Overlap Brace (BOB) Stock Size 15°

The Original Boston Overlap Brace (BOB) is an anterior opening LSO thermoplastic brace specifically designed for patients suffering from chronic or temporary low back pain.

Recommended Range Of Application

Mechanical back pain, Discongenic low back pain, Compression fractures, Sports injuries, Spondylolysis, Spondylolisthesis, Osteoporosis, Spinal Arthrosis, Spinal Stenosis, Pre and Post-operative care.

Measurements

Available in 20 stock sizes with either 0°, 15° or 30° Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree
151	Stock Size Module, Unlined	1	71	51	15°
151D	Stock Size Module, Unlined	1D	71	45	15°
152	Stock Size Module, Unlined	2	76	61	15°
152D	Stock Size Module, Unlined	2D	76	51	15°
153	Stock Size Module, Unlined	3	81	71	15°
153D	Stock Size Module, Unlined	3D	81	66	15°
154	Stock Size Module, Unlined	4	86.5	76	15°
154D	Stock Size Module, Unlined	4D	86.5	58.5	15°
155	Stock Size Module, Unlined	5	91.5	81	15°
155D	Stock Size Module, Unlined	5D	91.5	66	15°
156	Stock Size Module, Unlined	6	96.5	86.5	15°
156D	Stock Size Module, Unlined	6D	96.5	74	15°
157	Stock Size Module, Unlined	7	101.5	91.5	15°
157D	Stock Size Module, Unlined	7D	101.5	81	15°
158	Stock Size Module, Unlined	8	109	96.5	15°
158D	Stock Size Module, Unlined	8D	109	84	15°
159	Stock Size Module, Unlined	9	114	100	15°
159D	Stock Size Module, Unlined	9D	114	87	15°
1510	Stock Size Module, Unlined	10	117	104	15°
1510D	Stock Size Module, Unlined	10D	117	96	15°
151L	Stock Size Module, Lined	1	71	51	15°
151DL	Stock Size Module, Lined	1D	71	45	15°
152L	Stock Size Module, Lined	2	76	45	15°
152DL	Stock Size Module, Lined	2D	76	51	15°
153L	Stock Size Module, Lined	3	81	71	15°
153DL	Stock Size Module, Lined	3D	81	66	15°
154L	Stock Size Module, Lined	4	86.5	76	15°
154DL	Stock Size Module, Lined	4D	86.5	58.5	15°
155L	Stock Size Module, Lined	5	91.5	81	15°
155DL	Stock Size Module, Lined	5D	91.5	66	15°
156L	Stock Size Module, Lined	6	96.5	86.5	15°
156DL	Stock Size Module, Lined	6D	96.5	74	15°
157L	Stock Size Module, Lined	7	101.5	91.5	15°
157DL	Stock Size Module, Lined	7D	101.5	81	15°
158L	Stock Size Module, Lined	8	109	96.5	15°
158DL	Stock Size Module, Lined	8D	109	84	15°
159L	Stock Size Module, Lined	9	114	100	15°
159DL	Stock Size Module, Lined	9D	114	87	15°
1510L	Stock Size Module, Lined	10	117	104	15°
1510DL	Stock Size Module, Lined	10D	117	96	15°