

Boston Overlap Brace (BOB) Stock Size 30°

The Original Boston Overlap Brace (BOB) is an anterior opening LSO thermoplastic brace specifically designed for patients suffering from chronic or temporary low back pain.

Recommended Range Of Application

Mechanical back pain, Discongenic low back pain, Compression fractures, Sports injuries, Spondylolysis, Spondylolisthesis, Osteoporosis, Spinal Arthrosis, Spinal Stenosis, Pre and Post-operative care.

Measurements

Available in 20 stock sizes with either 0°, 15° or 30° Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree
301	Stock Size Module, Unlined	1	71	51	30°
301D	Stock Size Module, Unlined	1D	71	45	30°
302	Stock Size Module, Unlined	2	76	61	30°
302D	Stock Size Module, Unlined	2D	76	51	30°
303	Stock Size Module, Unlined	3	81	71	30°
303D	Stock Size Module, Unlined	3D	81	66	30°
304	Stock Size Module, Unlined	4	86.5	76	30°
304D	Stock Size Module, Unlined	4D	86.5	58.5	30°
305	Stock Size Module, Unlined	5	91.5	81	30°
305D	Stock Size Module, Unlined	5D	91.5	66	30°
306	Stock Size Module, Unlined	6	96.5	86.5	30°
306D	Stock Size Module, Unlined	6D	96.5	74	30°
307	Stock Size Module, Unlined	7	101.5	91.5	30°
307D	Stock Size Module, Unlined	7D	101.5	81	30°
308	Stock Size Module, Unlined	8	109	96.5	30°
308D	Stock Size Module, Unlined	8D	109	84	30°
309	Stock Size Module, Unlined	9	114	100	30°
309D	Stock Size Module, Unlined	9D	114	87	30°
310	Stock Size Module, Unlined	10	117	104	30°
310D	Stock Size Module, Unlined	10D	117	96	30°
301L	Stock Size Module, Lined	1	71	51	30°
301DL	Stock Size Module, Lined	1D	71	45	30°
302L	Stock Size Module, Lined	2	76	45	30°
302DL	Stock Size Module, Lined	2D	76	51	30°
303L	Stock Size Module, Lined	3	81	71	30°
303DL	Stock Size Module, Lined	3D	81	66	30°
304L	Stock Size Module, Lined	4	86.5	76	30°
304DL	Stock Size Module, Lined	4D	86.5	58.5	30°
305L	Stock Size Module, Lined	5	91.5	81	30°
305DL	Stock Size Module, Lined	5D	91.5	66	30°
306L	Stock Size Module, Lined	6	96.5	86.5	30°
306DL	Stock Size Module, Lined	6D	96.5	74	30°
307L	Stock Size Module, Lined	7	101.5	91.5	30°
307DL	Stock Size Module, Lined	7D	101.5	81	30°
308L	Stock Size Module, Lined	8	109	96.5	30°
308DL	Stock Size Module, Lined	8D	109	84	30°
309L	Stock Size Module, Lined	9	114	100	30°
309DL	Stock Size Module, Lined	9D	114	87	30°
3010L	Stock Size Module, Lined	10	117	104	30°
3010DL	Stock Size Module, Lined	10D	117	96	30°