

# Boston Body Jacket

The Boston Body Jacket is designed to immobilize spinal patients from neck to hips. This TLSO / LSO promotes rapid fusion for shorter hospital stays, and is therefore the ideal post-operative brace for your patients. Manufactured from low density polyethylene with optional hypoallergenic foam liner.

## Recommended Range Of Application

Post-operative, Compression Fractures, Spinal Immobilization

## Measurements

12 Bi-valve stock sizes available with 15° Lordosis only. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Xyphoid Can be made to measure or cast. Please use our Universal Measurement Form for ordering customised braces which can be found in our download centre.

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree
BJF01	Female Stock Size Module, Unlined, Bi-Valve	X-Small	88-99cm	69-79cm	15°
BJF02	Female Stock Size Module, Unlined, Bi-Valve	Small	94-104cm	74-84cm	15°
BJF03	Female Stock Size Module, Unlined, Bi-Valve	Medium	99-109cm	79-89cm	15°
BJF04	Female Stock Size Module, Unlined, Bi-Valve	Large	104-114cm	89-99cm	15°
BJF05	Female Stock Size Module, Unlined, Bi-Valve	X-Large	109-119cm	91-102cm	15°
BJF06	Female Stock Size Module, Unlined, Bi-Valve	XX-Large	119-130cm	112-122cm	15°
BJF01-L	Female Stock Size Module, Lined, Bi-Valve	X-Small	88-99cm	69-79cm	15°
BJF02-L	Female Stock Size Module, Lined, Bi-Valve	Small	94-104cm	74-84cm	15°
BJF03-L	Female Stock Size Module, Lined, Bi-Valve	Medium	99-109cm	79-89cm	15°
BJF04-L	Female Stock Size Module, Lined, Bi-Valve	Large	104-114cm	89-99cm	15°
BJF05-L	Female Stock Size Module, Lined, Bi-Valve	X-large	109-119cm	91-102cm	15°
BJF06-L	Female Stock Size Module, Lined, Bi-Valve	XX-Large	119-130cm	112-122cm	15°
BJM01	Male Stock Size Module, Unlined, Bi-Valve	X-Small	89-99cm	76-86cm	15°
BJM02	Male Stock Size Module, Unlined, Bi-Valve	Small	94-104cm	81-91cm	15°
BJM03	Male Stock Size Module, Unlined, Bi-Valve	Medium	99-109cm	86-96cm	15°
BJM04	Male Stock Size Module, Unlined, Bi-Valve	Large	104-114cm	96-107cm	15°
BJM05	Male Stock Size Module, Unlined, Bi-Valve	X-Large	109-119cm	99-109cm	15°
BJM06	Male Stock Size Module, Unlined, Bi-Valve	XX-Large	119-130cm	119-130cm	15°
BJM01-L	Male Stock Size Module, Lined, Bi-Valve	X-Small	89-99cm	76-86cm	15°
BJM02-L	Male Stock Size Module, Lined, Bi-Valve	Small	94-104cm	81-91cm	15°
BJM03-L	Male Stock Size Module, Lined, Bi-Valve	Medium	99-109cm	86-96cm	15°
BJM04-L	Male Stock Size Module, Lined, Bi-Valve	Large	104-114cm	96-107cm	15°
BJM05-L	Male Stock Size Module, Lined, Bi-Valve	X-Large	109-119cm	99-109cm	15°
BJM06-L	Male Stock Size Module, Lined, Bi-Valve	XX-Large	119-130cm	119-130cm	15°
BBJ011	Made to Measure, Lined, Ant / Post / Lat / TLSO / LSO	Customised	-	-	Customised
BBJ00	Made to Measure Module, Lined, Bi-Valve / TLSO / LSO	Customised	-	-	Customised
BBJ001	Made to Measure Module, Unlined, Bi-Valve / TLSO / LSO	Customised	-	-	Customised
BBJ033	Made to Cast Module, Lined, Ant / Post / Lat / TLSO / LSO	Customised	-	-	Customised
BBJ02	Made to Cast Module, Lined, Bi-Valve / TLSO / LSO	Customised	-	-	Customised
BBJ002	Made to Cast, Unlined, Bi-Valve / TLSO / LSO	Customised	-	-	Customised
B-FINISH	Finishing to First Fitting	Customised	-	-	Customised