

Boston Soft Overlap Brace (Soft BOB)

The Soft Boston Overlap Brace (Soft BOB) is an anterior opening brace made from 3/16 inch low density thermofoam and a 1/8 inch high density thermofoam. This is a versatile, easy to apply medium duty method for treating most indications where back pain is evident.

Recommended Range Of Application

Relief of chronic or temporary severe back pain; treatment and prevention of sports injuries and treating specific spinal indications.

Measurements

Available in 20 stock sizes with either 0 or 15 degree Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest. Can also be manufactured to cast or measurements. Please use our universal measurement form for customised brace ordering.

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree	Measurement 3
SBB-001	Stock size module	1	71	51	0°	66
SBB-001D	Stock size module	1D	71	45	0°	58.5
SBB-002	Stock size module	2	76	45	0°	58.5
SBB-002D	Stock size module	2D	76	51	0°	63.5
SBB-003	Stock size module	3	81	71	0°	76
SBB-003D	Stock size module	3D	81	66	0°	74
SBB-004	Stock size module	4	86.5	76	0°	84
SBB-004D	Stock size module	4D	86.5	58.5	0°	69
SBB-005	Stock size module	5	91.5	81	0°	89
SBB-005D	Stock size module	5D	91.5	66	0°	76
SBB-006	Stock size module	6	96.5	86.5	0°	94
SBB-006D	Stock size module	6D	96.5	74	0°	84
SBB-007	Stock size module	7	101.5	91.5	0°	99
SBB-007D	Stock size module	7D	101.5	81	0°	91.5
SBB-008	Stock size module	8	109	96	0°	109
SBB-008D	Stock size module	8D	109	84	0°	91.5
SBB-009	Stock size module	9	114	100	0°	114
SBB-009D	Stock size module	9D	114	87	0°	94
SBB-0010	Stock size module	10	117	104	0°	117
SBB-0010D	Stock size module	10D	117	96	0°	101
SBB-151	Stock size module	1	71	51	15°	66
SBB-151D	Stock size module	1D	71	45	15°	58.5
SBB-152	Stock size module	2	76	45	15°	58.5
SBB-152D	Stock size module	2D	76	51	15°	63.5
SBB-153	Stock size module	3	81	71	15°	76
SBB-153D	Stock size module	3D	81	66	15°	74
SBB-154	Stock size module	4	86.5	76	15°	84
SBB-154D	Stock size module	4D	86.5	58.5	15°	69
SBB-155	Stock size module	5	91.5	81	15°	89
SBB-155D	Stock size module	5D	91.5	66	15°	76
SBB-156	Stock size module	6	96.5	86.5	15°	94
SBB-156D	Stock size module	6D	96.5	74	15°	84
SBB-157	Stock size module	7	101.5	91.5	15°	99
SBB-157D	Stock size module	7D	101.5	81	15°	91.5
SBB-158	Stock size module	8	109	96	15°	109
SBB-158D	Stock size module	8D	109	84	15°	91.5
SBB-159	Stock size module	9	114	100	15°	114
SBB-159D	Stock size module	9D	114	87	15°	94
SBB-1510	Stock size module	10	117	104	15°	117
SBB-1510D	Stock size module	10D	117	96	15°	101
SBOB04	Made to Measure Module, Lined	Customised	-	-	Customised	-