

# Boston Soft Body Jacket Customised

The original Soft Boston Body Jacket (BSBJ) is a spinal orthosis designed to make control more comfortable for patients who cannot tolerate rigid compression. The BSBJ remains flexible while providing support with strategically placed co-polymer Stays (reinforcements) which increase diaphragm expansion to improve respiration. Choose from Removable or Permanent Stays. Alternatively, the product can be manufactured with either an internal or external polypropylene frame, where slightly more rigidity is required.

## Recommended Range Of Application

Neuromuscular Scoliosis, Seating and Positioning (posture control), Myelomeningocele, Cerebral Palsy, Muscular Dystrophy. Especially beneficial for children and non-surgical candidates.

## Measurements

Can be manufactured to cast or measurement with posterior, anterior or bi-valve opening. Please use our Universal Measurement Form for ordering customised braces which can be found in our download centre.

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.



Item No.	Model	Measurement	Lordosis degree
BSBJ02	Made to Measure Module, Removable Stays	Customised	Customised
BSBJ03	Made to Measure Module, Permanent Stays	Customised	Customised
BSBJ04	Made to Cast Module, Removable Stays	Customised	Customised
BSBJ05	Made to Cast Module, Permanent Stays	Customised	Customised
BSBJ06	Made to Cast Module, Internal/External Frame	Customised	Customised
BSBJ07	Made to Measure Module, Internal/External Frame	Customised	Customised
BSBJ/BI	Bi-valve Opening on Soft Body Jacket	-	-
BSBJ/COL	Coloured Foam - Blue/Orange/Purple/Pink	-	-
BSBJ/FIN	Finishing to First Fitting	-	-
BSBJ/HS	Head Support / Outtrigger	Customised	-