

Boston Soft Body Jacket Stock Size 0° Removable Stays



The original Soft Boston Body Jacket (BSBJ) is a spinal orthosis designed to make control more comfortable for patients who cannot tolerate rigid compression. The BSBJ remains flexible while providing support with strategically placed co-polymer Stays (reinforcements) which increase diaphragm expansion to improve respiration. Choose from Removable or Permanent Stays.

Recommended Range Of Application

Neuromuscular Scoliosis, Seating and Positioning (posture control), Myelomeningocele, Cerebral Palsy, Muscular Dystrophy. Especially beneficial for children and non-surgical candidates.

Measurements

Available in 30 stock sizes which can be manufactured with either an anterior or posterior opening. Supplied in either 0° or 15° Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.

Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree	Measurement 3
SBB-00J1R	Stock Size Module, Removable Stays	J1	54	40	0°	49
SBB-00J2R	Stock Size Module, Removable Stays	J2	58	42	0°	51
SBB-00J3R	Stock Size Module, Removable Stays	J3	62	45	0°	54
SBB-00J4R	Stock Size Module, Removable Stays	J4	64	47	0°	56
SBB-00Y1AR	Stock Size Module, Removable Stays	Y1A	66	51	0°	55
SBB-00Y1R	Stock Size Module, Removable Stays	Y1	67	46	0°	56
SBB-00Y2R	Stock Size Module, Removable Stays	Y2	71	47	0°	60
SBB-00Y2AR	Stock Size Module, Removable Stays	Y2A	72	55	0°	62
SBB-00Y3R	Stock Size Module, Removable Stays	Y3	75	48	0°	63
SBB-00Y3AR	Stock Size Module, Removable Stays	Y3A	78	57	0°	69
SBB-00Y4R	Stock Size Module, Removable Stays	Y4	79	51	0°	64
SBB-00Y4FR	Stock Size Module, Removable Stays	Y4F	80	57	0°	66
SBB-00Y5R	Stock Size Module, Removable Stays	Y5	82	52	0°	68
SBB-00Y5FR	Stock Size Module, Removable Stays	Y5F	82	58	0°	68
SBB-00Y6R	Stock Size Module, Removable Stays	Y6	84	56	0°	73
SBB-00Y6AR	Stock Size Module, Removable Stays	Y6A	84	63	0°	78
SBB-00Y7R	Stock Size Module, Removable Stays	Y7	86	54	0°	69
SBB-00Y7FR	Stock Size Module, Removable Stays	Y7F	86	60	0°	69
SBB-00Y7AR	Stock Size Module, Removable Stays	Y7A	86	69	0°	73
SBB-00Y8R	Stock Size Module, Removable Stays	Y8	88	58	0°	75
SBB-00Y8FR	Stock Size Module, Removable Stays	Y8F	88	64	0°	75
SBB-00Y8AR	Stock Size Module, Removable Stays	Y8A	90	72	0°	83
SBB-00Y9R	Stock Size Module, Removable Stays	Y9	92	60	0°	72
SBB-00Y9AR	Stock Size Module, Removable Stays	Y9A	92	63	0°	69
SBB-00Y9FR	Stock Size Module, Removable Stays	Y9F	92	68	0°	78
SBB-00A1R	Stock Size Module, Removable Stays	A1	94	62	0°	76
SBB-00A2R	Stock Size Module, Removable Stays	A2	94	70	0°	84
SBB-00A2AR	Stock Size Module, Removable Stays	A2A	96	66	0°	75
SBB-00A3R	Stock Size Module, Removable Stays	A3	97	68	0°	80