

Boston Soft Body Jacket Stock Size 0° Permanent Stays



The original Soft Boston Body Jacket (BSBJ) is a spinal orthosis designed to make control more comfortable for patients who cannot tolerate rigid compression. The BSBJ remains flexible while providing support with strategically placed co-polymer Stays (reinforcements) which increase diaphragm expansion to improve respiration. Choose from Removable or Permanent Stays.

Recommended Range Of Application

Neuromuscular Scoliosis, Seating and Positioning (posture control), Myelomeningocele, Cerebral Palsy, Muscular Dystrophy. Especially beneficial for children and non-surgical candidates.

Measurements

Available in 30 stock sizes which can be manufactured with either an anterior or posterior opening. Supplied in either 0° or 15° Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest.

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.

Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree	Measurement 3
SBB-00J1P	Stock Size Module, Permanent Stays	J1	54	40	0°	49
SBB-00J2P	Stock Size Module, Permanent Stays	J2	58	42	0°	51
SBB-00J3P	Stock Size Module, Permanent Stays	J3	62	45	0°	54
SBB-00J4P	Stock Size Module, Permanent Stays	J4	64	47	0°	56
SBB-00Y1AP	Stock Size Module, Permanent Stays	Y1A	66	51	0°	55
SBB-00Y1P	Stock Size Module, Permanent Stays	Y1	67	46	0°	56
SBB-00Y2P	Stock Size Module, Permanent Stays	Y2	71	47	0°	60
SBB-00Y2AP	Stock Size Module, Permanent Stays	Y2A	72	55	0°	62
SBB-00Y3P	Stock Size Module, Permanent Stays	Y3	75	48	0°	63
SBB-00Y3AP	Stock Size Module, Permanent Stays	Y3A	78	57	0°	69
SBB-00Y4P	Stock Size Module, Permanent Stays	Y4	79	51	0°	64
SBB-00Y4FP	Stock Size Module, Permanent Stays	Y4F	80	57	0°	66
SBB-00Y5P	Stock Size Module, Permanent Stays	Y5	82	52	0°	68
SBB-00Y5FP	Stock Size Module, Permanent Stays	Y5F	82	58	0°	68
SBB-00Y6P	Stock Size Module, Permanent Stays	Y6	84	56	0°	73
SBB-00Y6AP	Stock Size Module, Permanent Stays	Y6A	84	63	0°	78
SBB-00Y7P	Stock Size Module, Permanent Stays	Y7	86	54	0°	69
SBB-00Y7FP	Stock Size Module, Permanent Stays	Y7F	86	60	0°	69
SBB-00Y7AP	Stock Size Module, Permanent Stays	Y7A	86	69	0°	73
SBB-00Y8P	Stock Size Module, Permanent Stays	Y8	88	58	0°	75
SBB-00Y8FP	Stock Size Module, Permanent Stays	Y8F	88	64	0°	75
SBB-00Y8AP	Stock Size Module, Permanent Stays	Y8A	90	72	0°	83
SBB-00Y9P	Stock Size Module, Permanent Stays	Y9	92	60	0°	72
SBB-00Y9AP	Stock Size Module, Permanent Stays	Y9A	92	63	0°	69
SBB-00Y9FP	Stock Size Module, Permanent Stays	Y9F	92	68	0°	78
SBB-00A1P	Stock Size Module, Permanent Stays	A1	94	62	0°	76
SBB-00A2P	Stock Size Module, Permanent Stays	A2	94	70	0°	84
SBB-00A2AP	Stock Size Module, Permanent Stays	A2A	96	66	0°	75
SBB-00A3P	Stock Size Module, Permanent Stays	A3	97	68	0°	80
SBB-00A4P	Stock Size Module, Permanent Stays	A4	97	76	0°	88