

Boston Soft Body Jacket Stock Size 15° Permanent Stays



The original Soft Boston Body Jacket (BSBJ) is a spinal orthosis designed to make control more comfortable for patients who cannot tolerate rigid compression. The BSBJ remains flexible while providing support with strategically placed co-polymer Stays (reinforcements) which increase diaphragm expansion to improve respiration. Choose from Removable or Permanent Stays.

Recommended Range Of Application

Neuromuscular Scoliosis, Seating and Positioning (posture control), Myelomeningocele, Cerebral Palsy, Muscular Dystrophy. Especially beneficial for children and non-surgical candidates.

Measurements

Available in 30 stock sizes which can be manufactured with either an anterior or posterior opening. Supplied in either 0° or 15° Lordosis. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest
This product can also be made to scan (our new Boston Scan-to-Fit offering).
Follow the protocol for scanning and the instruction video for the correct procedure.
The information can be found under Support documents and Videos below.

Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree	Measurement 3
SBB-15J1P	Stock Size Module, Permanent Stays	J1	54	40	15°	49
SBB-15J2P	Stock Size Module, Permanent Stays	J2	58	42	15°	51
SBB-15J3P	Stock Size Module, Permanent Stays	J3	62	45	15°	54
SBB-15J4P	Stock Size Module, Permanent Stays	J4	64	47	15°	56
SBB-15Y1AP	Stock Size Module, Permanent Stays	Y1A	66	51	15°	55
SBB-15Y1P	Stock Size Module, Permanent Stays	Y1	67	46	15°	56
SBB-15Y2P	Stock Size Module, Permanent Stays	Y2	71	47	15°	60
SBB-15Y2AP	Stock Size Module, Permanent Stays	Y2A	72	55	15°	62
SBB-15Y3P	Stock Size Module, Permanent Stays	Y3	75	48	15°	63
SBB-15Y3AP	Stock Size Module, Permanent Stays	Y3A	78	57	15°	69
SBB-15Y4P	Stock Size Module, Permanent Stays	Y4	79	51	15°	64
SBB-15Y4FP	Stock Size Module, Permanent Stays	Y4F	80	57	15°	66
SBB-15Y5P	Stock Size Module, Permanent Stays	Y5	82	52	15°	68
SBB-15Y5FP	Stock Size Module, Permanent Stays	Y5F	82	58	15°	68
SBB-15Y6P	Stock Size Module, Permanent Stays	Y6	84	56	15°	73
SBB-15Y6AP	Stock Size Module, Permanent Stays	Y6A	84	63	15°	78
SBB-15Y7P	Stock Size Module, Permanent Stays	Y7	86	54	15°	69
SBB-15Y7FP	Stock Size Module, Permanent Stays	Y7F	86	60	15°	69
SBB-15Y7AP	Stock Size Module, Permanent Stays	Y7A	86	69	15°	73
SBB-15Y8P	Stock Size Module, Permanent Stays	Y8	88	58	15°	75
SBB-15Y8FP	Stock Size Module, Permanent Stays	Y8F	88	64	15°	75
SBB-15Y8AP	Stock Size Module, Permanent Stays	Y8A	90	72	15°	83
SBB-15Y9P	Stock Size Module, Permanent Stays	Y9	92	60	15°	72
SBB-15Y9AP	Stock Size Module, Permanent Stays	Y9A	92	63	15°	69
SBB-15Y9FP	Stock Size Module, Permanent Stays	Y9F	92	68	15°	78
SBB-15A1P	Stock Size Module, Permanent Stays	A1	94	62	15°	76
SBB-15A2P	Stock Size Module, Permanent Stays	A2	94	70	15°	84
SBB-15A2AP	Stock Size Module, Permanent Stays	A2A	96	66	15°	75
SBB-15A3P	Stock Size Module, Permanent Stays	A3	97	68	15°	80
SBB-15A4P	Stock Size Module, Permanent Stays	A4	97	76	15°	88