

Boston Lite™

Boston Lite™ is a new concept based on the Boston Soft Brace. New techniques have made it possible to create a brace that is both thinner and more ventilated, but yet with the same proven stability. The new brace has new benefits and functions. By merging the two materials and using certain shapes, it is easier to accommodate the brace to the patients need for ventilation or feeding tubes. The lining is a separate piece made from a 3D-spacer material. This material provides comfort by absorbing moisture and increases ventilation to keep the skin cooler and drier. The spacer also works as a pressure divider towards the body. The liner is used instead of a traditional t-shirt. To even further improve the cosmetic design, the new profile buckles have been developed that contributes to a more low profile design.



Recommended Range Of Application

Neuromuscular Scoliosis, Trunk control and Post-operative support.

Measurements

Available in 12 stock sizes. Can be manufactured to cast or measurement. Please use our Universal Measurement Form in our download centre for ordering custom braces. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.

Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree	Measurement 3
SBBL-15J1	Stock Size Module	J1	54	40	15°	49
SBBL-15J2	Stock Size Module	J2	58	42	15°	51
SBBL-15J3	Stock Size Module	J3	62	45	15°	54
SBBL-15J4	Stock Size Module	J4	64	47	15°	56
SBBL-15Y1	Stock Size Module	Y1	67	46	15°	56
SBBL-15Y1A	Stock Size Module	Y1A	66	51	15°	55
SBBL-15Y2	Stock Size Module	Y2	71	47	15°	60
SBBL-15Y2A	Stock Size Module	Y2A	72	55	15°	62
SBBL-15Y3	Stock Size Module	Y3	75	48	15°	63
SBBL-15Y3A	Stock Size Module	Y3A	78	57	15°	69
SBBL-15Y4	Stock Size Module	Y4	79	51	15°	64
SBBL-15Y4F	Stock Size Module	Y4F	80	57	15°	66
BSBJL02	Made to Measure Module	Customised	-	-	-	-
BSBJL03	Made to Cast Module	Customised	-	-	-	-
BSBJL/FIN	Finishing to First Fitting	Customised	-	-	-	-
SBBL-L1	Liner	1 (15J1-15J2)	-	-	-	-
SBBL-L2	Liner	2 (15J3-15J4)	-	-	-	-
SBBL-L3	Liner	3 (15Y1-15Y2)	-	-	-	-
SBBL-L4	Liner	4 (15Y1A-15Y2A)	-	-	-	-
SBBL-L5	Liner	5 (15Y3-15Y4)	-	-	-	-
SBBL-L6	Liner	6 (15Y3A-15Y4F)	-	-	-	-
SBBL-S	Single Strap & Buckle	One Size	-	-	-	-
SBBL-S/SET	Set of 3 Straps & Buckles	One Size	-	-	-	-
SBBL-B	Buckle Cover	One Size	-	-	-	-