

Boston Transitional Belt (BTB)



The Boston Transitional Belt (BTB) is designed for use during physical activity and sporting events. It is constructed with an inner copolymer frame sandwiched between a layer of 3/16 inch thermofoam on the inside and 1/8 inch high density thermofoam on the outside.

Recommended Range Of Application

The patient typically wears the BTB for a 12 month period following successful completion of treatment with a Boston Overlap Brace for Spondylolisthesis, Spondylolysis and mechanical low back pain.

Measurements

Available in 20 stock sizes. Supplied already finished. Can also be made to measurements. Please use our Universal Measurement Form for ordering customised braces which can be found in our download centre. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.

Item No.	Model	Size	Measurement 1	Measurement 2	Measurement 3
BTB-1	Stock Size Module	1	71	51	66
BTB-2	Stock Size Module	2	71	45	58.5
BTB-3	Stock Size Module	3	76	45	58.5
BTB-4	Stock Size Module	4	76	51	63.5
BTB-5	Stock Size Module	5	81	71	76
BTB-6	Stock Size Module	6	81	66	74
BTB-7	Stock Size Module	7	86.5	76	84
BTB-8	Stock Size Module	8	86.5	58.5	69
BTB-9	Stock Size Module	9	91.5	81	89
BTB-10	Stock Size Module	10	91.5	66	76
BTB-11	Stock Size Module	11	96.5	86.5	94
BTB-12	Stock Size Module	12	96.5	74	84
BTB-13	Stock Size Module	13	101.5	91.5	99
BTB-14	Stock Size Module	14	101.5	81	91.5
BTB-15	Stock Size Module	15	109	99.6	109
BTB-16	Stock Size Module	16	109	84	91.5
BTB-17	Stock Size Module	17	114	100	114
BTB-18	Stock Size Module	18	114	87	94
BTB-19	Stock Size Module	19	117	104	117
BTB-20	Stock Size Module	20	117	96	101
BTB-MTM	Made to Measure Module	Customised	-	-	-