

# Boston Combi Soft (BCS)

The Boston Combi Soft brace is an anterior opening TLSO. The semi-rigid structure combines stability with the inner copolymer frame, outer high density thermofoam and comfort of the soft inner thermofoam. The brace is designed for patients that would require a more traditional rigid body jacket, but are unable to tolerate the hard plastic type brace.

## Recommended Range Of Application

Osteoporosis, Spondylolysis / Spondylolisthesis, Compression Fractures and Post-Op Surgery.

## Measurements

Available in 10 stock sizes with 15° Lordosis. Can be manufactured to cast or measurement. Please use our Universal Measurement Form for ordering customised braces which can be found in our download centre. Measurements in centimetres. Measurement 1: Hip; Measurement 2: Waist; Measurement 3: Chest

This product can also be made to scan (our new Boston Scan-to-Fit offering). Follow the protocol for scanning and the instruction video for the correct procedure. The information can be found under Support documents and Videos below.



Item No.	Model	Size	Measurement 1	Measurement 2	Lordosis degree	Measurement 3
BCS-1	Stock Size Module	1	71	51	15°	66
BCS-2	Stock Size Module	2	76	45	15°	58.5
BCS-3	Stock Size Module	3	81	71	15°	76
BCS-4	Stock Size Module	4	86.5	76	15°	84
BCS-5	Stock Size Module	5	91.5	81	15°	89
BCS-6	Stock Size Module	6	96.5	86.5	15°	94
BCS-7	Stock Size Module	7	101.5	91.5	15°	99
BCS-8	Stock Size Module	8	109	96	15°	109
BCS-9	Stock Size Module	9	114	100	15°	114
BCS-10	Stock Size Module	10	117	104	15°	117
BCS-MTC	Made to Cast Module	Customised	-	-	Customised	-
BCS-MTM	Made to Measure Module	Customised	-	-	Customised	-