

Kiddie Composite Footplate Flat

A flat composite footplate that limits the range-of motion (ROM) in the meta-tarso-phalangeal joints (mtp-joints).

Recommended Range Of Application

The composite foot plate is intended to limit the range-of motion (ROM) and may be beneficial for Idiopathic toe walkers (ITW) who have difficulties with accepting ankle-foot-orthosis (AFO). Must be covered by an insole

Contraindications

Open ulcers on the foot. Non-stable fractures

Measurements

Length, Width, Size;



Item No.	Size	L/R	Length	Width	Thickness
280391011	Small	Left	175 mm	60 mm	1.6 mm
280391012	Medium	Left	190 mm	66 mm	1.6 mm
280391013	Large	Left	205 mm	72 mm	1.6 mm
280392011	Small	Right	175 mm	60 mm	1.6 mm
280392012	Medium	Right	190 mm	66 mm	1.6 mm
280392013	Large	Right	205 mm	72 mm	1.6 mm