

# SELECTION® Open Wrist Children

Selection Open Wrist Children is a brace that is simple to put on and take off. It can be opened totally flat to make it easier to put on. The integrated “wrap around” straps gives increased stability around the wrist. It has two thumb straps. One around the thumb and another on the dorsal side of the hand. This keeps the thumb in a correct position so that it does not drop towards the palm of the hand. The orthosis is intended to immobilize the CMC joint and wrist in order to reduce pain, to guide the wrist into a neutral position and to support a wrist after trauma/post-operatively if needed.

## Recommended Range Of Application

Rheumatoid Arthritis, CP.

### Measurements

- Measure 1: Wrist circumference
- Measure 2: Palm height
- Measure 3: Total length on the volar side

SELECTION®



Item No.	Model	Size	L/R	Color	Length	Measurement 1	Measurement 2
353061011	Children	Small	Left	Black	14,5 cm	11,5 - 14 cm	7 - 8 cm
353061012	Children	Medium	Left	Black	17 cm	12 - 15,5 cm	8 - 9 cm
353061013	Children	Large	Left	Black	17,5 cm	14 - 16,5 cm	8 - 9 cm
353062011	Children	Small	Right	Black	14,5 cm	11,5 - 14 cm	7 - 8 cm
353062012	Children	Medium	Right	Black	17 cm	12 - 15,5 cm	8 - 9 cm
353062013	Children	Large	Right	Black	17,5 cm	14 - 16,5 cm	8 - 9 cm