

# ToeOFF® Flow 2 ½

Using our experience and development in manufacturing technology, together with your feedback, we have developed a new Allard AFO — ToeOFF® Flow 2 ½.

ToeOFF® Flow 2 ½ has a unique formula with new features and benefits which include:  
½ heel height of 7 mm in all sizes, more ROM in A-P plane and smoother surface underneath the footplate. The low profile tibia plate is more dynamic in the A-P plane, and it has a better fit in today's lower heeled shoes.

ToeOFF Flow 2½ has all Microfix pre-applied which makes it quicker to attach a SoftKIT™ 2.0, ComfortKIT™ 2.0 or CoverKIT™2.0. Closure straps are the same as for ToeOFF®2.0

## Recommended Range Of Application

Foot drop combined with moderate to no spasticity. Limb proprioception deficit and mild proximal deficit. For further information on uses and contraindications see AFO Professional Instruction.

## Contraindications

Severe edema, leg ulcers, moderate to severe spasticity, medial-lateral instability.

ToeOFF® Flow 2 ½



Item No.	Model	Size	L/R	Footplate Length	Height
289601010	ToeOFF® Flow2½	X-Small	Left	21 cm	36 cm
289601011	ToeOFF® Flow2½	Small	Left	23 cm	38 cm
289601012	ToeOFF® Flow2½	Medium	Left	24,5 cm	40,5 cm
289601013	ToeOFF® Flow2½	Large	Left	27 cm	43 cm
289601014	ToeOFF® Flow2½	X-Large	Left	28,5 cm	45 cm
289602010	ToeOFF® Flow2½	X-Small	Right	21 cm	36 cm
289602011	ToeOFF® Flow2½	Small	Right	23 cm	38 cm
289602012	ToeOFF® Flow2½	Medium	Right	24,5 cm	40,5 cm
289602013	ToeOFF® Flow2½	Large	Right	27 cm	43 cm
289602014	ToeOFF® Flow2½	X-Large	Right	28,5 cm	45 cm