

ToeOFF® Flow 2 ½

Using our experience and development in manufacturing technology, together with your feedback, we have developed a new Allard AFO – ToeOFF® Flow 2 ½.

ToeOFF® Flow 2 ½ has a unique formula with new features and benefits which include:
½ heel height of 7 mm in all sizes, more ROM in A-P plane and smoother surface underneath the footplate. The low profile tibia plate is more dynamic in the A-P plane, and it has a better fit in today’s lower heeled shoes.

ToeOFF Flow 2½ has all Microfix pre-applied which makes it quicker to attach a SoftKIT™ 2.0, ComfortKIT™ 2.0 or CoverKIT™2.0. Closure straps are the same as for ToeOFF®2.0

Recommended Range Of Application

Foot drop combined with moderate to no spasticity. Limb proprioception deficit and mild proximal deficit. For further information on uses and contraindications see AFO Professional Instruction.

Contraindications

Severe edema, leg ulcers, moderate to severe spasticity, medial-lateral instability.

ToeOFF® Flow 2 ½



Item No.	Model	Size	L/R	Footplate Length	Height
289601010	ToeOFF® Flow2½	X-Small	Left	21 cm	36 cm
289601011	ToeOFF® Flow2½	Small	Left	23 cm	38 cm
289601012	ToeOFF® Flow2½	Medium	Left	24,5 cm	40,5 cm
289601013	ToeOFF® Flow2½	Large	Left	27 cm	43 cm
289601014	ToeOFF® Flow2½	X-Large	Left	28,5 cm	45 cm
289602010	ToeOFF® Flow2½	X-Small	Right	21 cm	36 cm
289602011	ToeOFF® Flow2½	Small	Right	23 cm	38 cm
289602012	ToeOFF® Flow2½	Medium	Right	24,5 cm	40,5 cm
289602013	ToeOFF® Flow2½	Large	Right	27 cm	43 cm
289602014	ToeOFF® Flow2½	X-Large	Right	28,5 cm	45 cm