

# HyMo Textile Orthosis

HyMo is a dynamic textile orthosis for people with hEDS (hypermobile Ehlers-Danlos syndrome) or hypermobility within HSD (Hypermobility Spectrum Disorders). The made-to-measure orthosis provides even and constant pressure against the body, which may help improve posture, increase joint stability, reduce pain, and enhance proprioception.

HyMo is made from a high-quality elastic material that adapts to the body's shape and movements. The fabric is produced using a patented manufacturing process that results in a strong yet lightweight material that is comfortable to wear and maintains its stretch over time, even with daily use.

The garments are made-to-measure with a focus on stabilising the joints and, when needed, limiting the range of motion — for example in the hip or shoulder joints. When required, the orthosis can be supplemented with additional reinforcements for increased stability.

HyMo is designed as a bodysuit and/or trousers with short or long sleeves and legs, and can be used by adults and children from the age of 10. The orthosis is used gradually up to 8–10 hours per day. The materials are certified according to OEKO TEX® Standard 100, class II.

## Recommended Range Of Application

hEDS – Ehlers-Danlos Syndrome that presents with hypermobility, HSD (Hypermobility Spectrum Disorders)

## Contraindications

Oedema, dislocated joints (without reduction), Vascular EDS (vEDS).



Item No.	Model	Color	Type
394000000	Body with short sleeves	Black	Female
394010000	Body with short sleeves	Black	Male
394020000	Body with long sleeves	Black	Female
394030000	Body with long sleeves	Black	Male
394040000	Trousers with short legs	Black	Female
394050000	Trousers with short legs	Black	Male
394060000	Trousers with long legs	Black	Female
394070000	Trousers with long legs	Black	Male
394080000	Body with 3/4 sleeves	Black	Female
394090000	Body with 3/4 sleeves	Black	Male
394100000	Trousers with 3/4 legs	Black	Female
394110000	Trousers with 3/4 legs	Black	Male