

Elements Body Trial Sleeve

The Trial sleeve has been designed to be used in the assessment process for Elements Body MTM upper limb garments. It enables the clinician to evaluate the user's response to the application of a dynamic textile garment in preparation for a made to measure upper limb garment. The Trial sleeve is only to be used for an assessment period of 2–4 weeks to assist the clinician in identifying the clinical effectiveness and tolerance prior to the provision of an Elements Body MTM Garment.

Measurements

- Measure 1: Circumference wrist cm
- Measure 2: Circumference elbow cm
- Measure 3: Circumference biceps cm



Item No.	Model	Size	Measurement 1	Measurement 2	Measurement 3
393440008	EB Trial Sleeve	XXX-Small	11-12 cm	16-18 cm	16-18 cm
393440009	EB Trial Sleeve	XX-Small	12-13 cm	18-20 cm	18-20 cm
393440010	EB Trial Sleeve	X-Small	13-14 cm	20-22,5 cm	20-22 cm
393440011	EB Trial Sleeve	Small	14-15 cm	22,5-25 cm	21-23 cm
393440012	EB Trial Sleeve	Medium	14-16 cm	24-27 cm	22-25 cm
393440013	EB Trial Sleeve	Large	16-18 cm	27-30 cm	25-28 cm
393440014	EB Trial Sleeve	X-Large	18-20 cm	30-33 cm	28-31 cm
393440015	EB Trial Sleeve	XX-Large	20-22 cm	33-37 cm	31-34 cm