

SoftKIT™ Short

Every ToeOFF® should have a soft interface between the orthosis and the tibial crest. SoftKIT is a padding system that includes two narrow neoprene pads to create a pressure-free channel to help relieve any potential discomfort from bony prominences. The kit includes pre-cut pressure-sensitive MikroFIX tabs to easily affix the pad to ToeOFF®. An extra pad is included, so one can be worn while the other is laundered.



Item No.	Model	Size	Color
283910011	SoftKIT™ Short	Small	Black
283910012	SoftKIT™ Short	Medium	Black
283910013	SoftKIT™ Short	Large	Black
283560011	MicroFIX	Small	Black
283560012	MicroFIX	Medium	Black
283560013	MicroFIX	Large	Black
283570000	MicroFIX Ears	One size	Black