

# Elements Body Sock

Elements Body Sock is a dynamic textile ankle foot orthosis, intended to give light and even compression. The below the knee sock is made to measure and can be combined with a variety of reinforcements and colour options. The various reinforcements are designed to assist correct positioning of the foot, as well as provide sensory and proprioceptive feedback. The sock has a zip on the anterior, and closed toes.

## Recommended Range Of Application

Cerebral Palsy: hypotonia, spasticity, ataxia

MS, Stroke & Acquired Brain Injury: spasticity, ataxia

Muscular Dystrophy, Acquired Spinal Cord Injury and Down's Syndrome: hypotonia

Ataxia: spastic ataxia

## Contraindications

Ulcers in foot and/or lower leg. Circulatory disorders.



Item No.	Model
----------	-------

393250000	Sock
-----------	------