

BlueROCKER® Short

BlueROCKER® Short is specifically developed to accommodate individuals with a shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area. Tibia padding should always be worn for patient comfort.

Both models have a heel height of 7 mm in all sizes.

Recommended Range Of Application

Foot drop in combination with weak plantar flexors, and in some cases limited knee control. This may also be used for bilateral problems and partial foot amputations, most proximal level is Chopart.

Contraindications

Should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits or severe spasticity.

BlueROCKER® Short



Item No.	Model	Size	L/R	Footplate Length	Height
289951011	BlueROCKER® Short	Small	Left	21,5 cm	29,5 cm
289951012	BlueROCKER® Short	Medium	Left	23,5 cm	32 cm
289951013	BlueROCKER® Short	Large	Left	25,5 cm	34 cm
289952011	BlueROCKER® Short	Small	Right	21,5 cm	29,5 cm
289952012	BlueROCKER® Short	Medium	Right	23,5 cm	32 cm
289952013	BlueROCKER® Short	Large	Right	25,5 cm	34 cm