

SELECTION® Ulnar Deviation

SELECTION® Ulnar Deviation Splint is intended to improve finger alignment and enhance hand function and functional grip. This hand splint is exclusively designed and produced using our unique three layer Selection fabric. With its low profile design, it meets the contours of the hand, without impeding function. The fabric can be trimmed without fraying, and therefore allows further customization if required. The easy fastening surface allows the hook to be placed anywhere on the splint to optimise the directional pull needed. The Ulnar Deviation splint supports the arches of the hand whilst providing re-alignment, stabilisation and can contribute to the management of pain. Should you need to provide wrist stabilisation for your patient, this splint can also be worn in combination with a Selection wrist splint.



Recommended Range Of Application

For patients with inflammatory and degenerative diseases, such as rheumatoid arthritis, that has resulted in ulnar deviation. Selection Ulnar Deviation can also be used for hypermobile joints and weak muscles to provide stability. Not intended for use with fixed contractures in the hand.

Measurements

Circumference MCP 2-5

Item No.	Size	Measurement	L/R	Color
352181010	X-Small	13-16 cm	Left	Black
352181011	Small	16-19 cm	Left	Black
352181012	Medium	19-22 cm	Left	Black
352181013	Large	22-25 cm	Left	Black
352182010	X-Small	13-16 cm	Right	Black
352182011	Small	16-19 cm	Right	Black
352182012	Medium	19-22 cm	Right	Black
352182013	Large	22-25 cm	Right	Black