

# KiddieFLOW™ Trial Brace

This product is permanently labelled “Trial Brace”. It is only to be used for patient assessment, to determine the correct product and/or size to achieve optimum gait for the patient. We recommend the first assessment to be done in parallel bars. Includes SoftKIT - Applied!

## Recommended Range Of Application

KiddieFLOW™ Trial Brace | Foot Drop, Children



Item No.	Model	Size	L/R	Footplate Length	Height
288941011	KiddieFLOW Trial Brace	Small	Left	16 cm	22 cm
288941012	KiddieFLOW Trial Brace	Medium	Left	18 cm	25,7 cm
288941013	KiddieFLOW Trial Brace	Large	Left	20 cm	29,5 cm
288941014	KiddieFLOW Trial Brace	X-Large	Left	21 cm	31,5 cm
288942011	KiddieFLOW Trial Brace	Small	Right	16 cm	22 cm
288942012	KiddieFLOW Trial Brace	Medium	Right	18 cm	25,7 cm
288942013	KiddieFLOW Trial Brace	Large	Right	20 cm	29,5 cm
288942014	KiddieFLOW Trial Brace	X-Large	Right	21 cm	31,5 cm
288951011	KiddieFLOW w/D-ring Youth Trial Brace	Small	Left	16 cm	22 cm
288951012	KiddieFLOW w/D-ring Youth Trial Brace	Medium	Left	18 cm	25,7 cm
288951013	KiddieFLOW w/D-ring Youth Trial Brace	Large	Left	20 cm	29,5 cm