

TRUalign™ OA

TRUalign™ OA is a time-tested single upright OActive™ brace design in a patient preferred wrap-style brace. The Q-hinge technology provides significant angular correction to relieve pressure and pain of knee osteoarthritis.

Features and Benefits

Increased patient compliance as it is:

- Easy to apply.
- Breathable, comfortable and conforms to the anatomy.
- Q-hinge angular adjustment delivers optimized 3-point leverage for maximum relief.

Measurements

Measure 1: 15 cm above patella
Measure 2: 15 cm below patella



Item No.	Size	L/R	Knee Measurement 1	Knee Measurement 2
550181011	Small / Medium	Left	36-47cm	28-36cm
550181012	Medium / Large	Left	47-65cm	36-48cm
550181013	Large / X-large	Left	65-81cm	48-64cm
550182011	Small / Medium	Right	36-47cm	28-36cm
550182012	Medium / Large	Right	47-65cm	36-48cm
550182013	Large / X-large	Right	65-81cm	48-64cm