

# X-Static Boston T

The X-Static Boston T has the same benefits as the Standard T and is also available in double axilla flap. The material content in the single axilla flap comprises 49% polyester, 46.5% Nylon, 5%. X-Static® and 2% Lycra®. X-Static® fibres help transport moisture away from the skin and prevent odour. It also has thermal and anti-static properties to keep you cool and comfortable. The double axilla flap comprises 95.6% Polyester/X-Static® blend and 4.4% Lycra®.



Item No.	Model	Size	Length	Width	Circumference
1BX1BVML	X-Static Boston T	Medium, Single Axilla Flap	533 mm	203 mm	61-81 cm
1BX1BVSL	X-Static Boston T	Small, Single Axilla Flap	483 mm	152 mm	40-61 cm
1BX1BVLL	X-Static Boston T	Large, Single Axilla Flap	711 mm	254 mm	81-101 cm
1BX1BVXX	X-Static Boston T	X-Large, Single Axilla Flap	711 mm	305 mm	101-122 cm
1XBBDML	X-Static Boston T	Medium, Double Axilla Flap	533 mm	203 mm	61-81 cm
1XBBDL	X-Static Boston T	Small, Double Axilla Flap	483 mm	152 mm	40-61 cm
1XBBDLL	X-Static Boston T	Large, Double Axilla Flap	711 mm	254 mm	81-101 cm
1XBBDXX	X-Static Boston T	X-Large, Double Axilla Flap	711 mm	305 mm	101-122 cm