

Elements Body - Trial Lumbar Wrap

The Trial Lumbar Wrap is an evaluation tool to be used in the assessment process for Elements Body. It enables the clinician to evaluate the patient's response to the application of a dynamic textile garment in preparation for a made-to-measure garment. The Trial Lumbar Wrap is intended for short-term use for a maximum of two weeks.



Item No.	Size	Color	Pelvic Circumference
393460009	XX-Small	Black	50-60 cm
393460010	X-Small	Black	60-75 cm
393460011	Small	Black	75-90 cm
393460012	Medium	Black	90-105 cm
393460013	Large	Black	105-125 cm