

Carbon Footplate

The carbon footplate is intended to protect the MTP joints and the forefoot by limiting the ROM, distributing the pressure in the forefoot, and facilitating rollover. The footplate can also be used to gradually increase movement after immobilizing treatment with a post-op shoe or walker. It can improve balance after toe amputation, dig IV, and/or V.

The footplate should be covered by an insole or custom foot orthotic.



Recommended Range Of Application

Hallux rigidus, Hallux Limitus, stress fractures and other injuries or arthritis conditions in the forefoot/midfoot, forefoot trauma, sprained toe.

Item No.	Size	L/R	Length	Width	Thickness
280351011	Small	Left	235mm	77mm	1,8 mm
280351013	Large	Left	275mm	90mm	1,8 mm
280351014	X-Large	Left	305mm	99mm	1,8 mm
280352011	Small	Right	235mm	77mm/	1,8 mm
280352012	Medium	Right	255mm	80mm	1,8 mm
280352013	Large	Right	275mm	90mm	1,8 mm
280352014	X-Large	Right	305mm	99mm	1,8 mm