

Surestep GO SMO

SureStep GO is a prefabricated SMO, designed to manage mild to severe pronation associated with hypotonia. This SMO uses circumferential compression to stabilize the foot and ankle. The open heel allows more of the foot to “feel the floor”, which aids in sensory input and posture. For maximum function, this SMO should always be worn with suitable footwear.

A SureStep SMO is highly recommended to be used together with our KiddieFLOW™, KiddieGAIT® or KiddieROCKER® to encourage heel-to-toe gait, especially for patients with low tone and pronation.

Recommended Range Of Application

Mild to severe pronation, hypotonia, triplanar instability in weight bearing, inability to stand independently, mild toe walking.

Contraindications

Hypertonia, spasticity, severe toe walking.



Item No.	Model	Size	L/R	Package	Shoe Size	Unit/Pack
287422007	SureStep GO SMO	7T	Right	Unit	UK 6 / EU 23	-
287422008	-	8T	Right	Unit	UK 7 / EU 24	-
287422009	SureStep GO SMO	9T	Right	Unit	UK 8.5 / EU 26	-
287432010	SureStep GO SMO	10Y	Right	Unit	UK 9 / EU 27	-
287432011	SureStep GO SMO	11Y	Right	Unit	UK 10 / EU 28	-
287432012	SureStep GO SMO	12Y	Right	Unit	UK 11 / EU 29	-
287422006	SureStep GO SMO	6T	Right	Unit	UK 5.5 / EU 22	-
287422005	-	-	-	-	-	-
287390000	-	-	-	-	-	-
287422004	SureStep GO SMO	4T	Right	Unit	UK 4 / EU 20	-
287422003	SureStep GO SMO	3T	Right	Unit	UK 2 / EU 19	-
287431012	SureStep GO SMO	12Y	Left	Unit	UK 11 / EU 29	-
287431011	SureStep GO SMO	11Y	Left	Unit	UK 10 / EU 28	-
287431010	SureStep GO SMO	10Y	Left	Unit	UK 9 / EU 27	-
287421009	SureStep GO SMO	9T	Left	Unit	UK 8.5 / EU 26	-
287421008	SureStep GO SMO	8T	Left	Unit	UK 7 / EU 24	-
287421007	SureStep GO SMO	7T	Left	Unit	UK 6 / EU 23	-
287421006	SureStep GO SMO	6T	Left	Unit	UK 5.5 / EU 22	-
287421005	SureStep GO SMO	5T	Left	Unit	UK 5 / EU 21	-
287421003	SureStep GO SMO	3T	Left	Unit	UK 2 / EU 19	-
287421004	SureStep GO SMO	4T	Left	Unit	UK 4 / EU 20	-
287380000	-	-	-	-	-	-
287400003	SureStep GO SMO	-	-	-	EU 19	1 pair
287400004	SureStep GO SMO	-	-	-	EU 20	1 pair
287400005	SureStep GO SMO	-	-	-	EU 21	1 pair
287400006	SureStep GO SMO	-	-	-	EU 22	1 pair
287400007	SureStep GO SMO	-	-	-	EU 23	1 pair
287400008	SureStep GO SMO	-	-	-	EU 24	1 pair
287400009	SureStep GO SMO	-	-	-	EU 26	1 pair
287410010	SureStep GO SMO	-	-	-	EU 27	1 pair
287410011	SureStep GO SMO	-	-	-	EU 28	1 pair
287410012	SureStep GO SMO	-	-	-	EU 29	1 pair