ToeOFF® Short Trial Brace



This product is permanently labelled “Trial Brace”. It is only to be used for patient assessment, to determine ToeOFF® product and/or size to achieve optimum gait for the patient. We recommend the first assessment to be done in parallel bars. Includes SoftKIT — Applied!

Recommended Range Of Application

This product is intended only for in-clinic trial assessment only.

Item No.	Model	Size	L/R	Footplate Length	Height
289911011	ToeOFF® Short Trial Brace	Small	Left	21,5 cm	29,5 cm
289911012	ToeOFF® Short Trial Brace	Medium	Left	23,5 cm	32 cm
289911013	ToeOFF® Short Trial Brace	Large	Left	25,5 cm	34 cm
289912011	ToeOFF® Short Trial Brace	Small	Right	21,5 cm	29,5 cm
289912012	ToeOFF® Short Trial Brace	Medium	Right	23,5 cm	32 cm
289912013	ToeOFF® Short Trial Brace	Large	Right	25,5 cm	34 cm